

# T.B.R.<sup>®</sup>

## Switch

NEW



### T.B.R.<sup>®</sup> Switch Abutments :

- Push away the biological space at the level of the implantary platform
- Protect the crestal bone level during 2 stages surgeries
- Keep the osseous volume
- Do not expose any more crestal bone to bacterial inflammation

### References :

	0°	15°	25°
Titanium screw-retained abutment for internal octagonal implant Ø 3.5 mm	O-MTS320	O-MTS321	O-MTS322
Titanium screw-retained abutment for internal octagonal implant Ø 4 mm	O-MTS420	O-MTS421	O-MTS422
Titanium screw-retained abutment for internal octagonal implant Ø 5 mm	O-MTS520	O-MTS521	O-MTS522



### **T.B.R.<sup>®</sup> GROUP**

24, impasse René Couzinet 31500 TOULOUSE / FRANCE  
 Tel: (+33) (0)5 62 16 71 00 / Fax: (+33) (0)5 61 80 84 02  
 e-mail: info@tbr-group.com

[www.tbr-group.com](http://www.tbr-group.com)



The smile - architects